



# A proximity focusing RICH with multiple refractive index aerogel radiator

#### Peter Križan University of Ljubljana and J. Stefan Institute

#### For Belle Aerogel RICH R&D group

October 20, 2005

9th ICATPP Conference, Como

Peter Križan, Ljubljana







Introduction, motivation and requirements Radiator with multiple refractive indices Beam test results Further radiator optimisation R&D issues



#### Accumulated data sample $\cong$ 500 M BB-pairs



improve K/ $\pi$  separation in the forward (high mom.) region for few-body decays of B's good K/ $\pi$  separation for b -> d $\gamma$ , b -> s $\gamma$ improve purity in fully reconstructed B decays low momentum (<1GeV/c) e/ $\mu$ / $\pi$  separation (B ->KII) keep high the efficiency for tagging kaons

9<sup>th</sup> ICATPP Conference, Como





Adachi, S. Fratina, T. Fukushima, A. Gorišek, T. Iijima, H. Kawai,
M. Konishi, S. Korpar, Y. Kozakai, P. Križan, T. Matsumoto,
S. Nishida, S. Ogawa, S. Ohtake, R. Pestotnik, S. Saitoh, T. Seki,
A. Stanovnik, T. Sumiyoshi, Y. Uchida, Y. Unno, S. Yamamoto

Chiba University, Japan Toho University, Funabashi, Japan High Energy Accelerator Research Organization (KEK), Japan J. Stefan Institute, Ljubljana, Slovenia University of Ljubljana, Slovenia Nagoya University, Nagoya, Japan Tokyo Metropolitan University



# Proximity focusing RICH in the forward region





K/π separation at 4 GeV/c  $\theta_c(\pi) \sim 308 \text{ mrad} (n = 1.05)$  $\theta_c(\pi) - \theta_c(K) \sim 23 \text{ mrad}$ 

 $d\theta_c$ (meas.) =  $\sigma_0 \sim 13$  mrad With 20mm thick aerogel and 6mm PMT pad size

 $\rightarrow$  6 $\sigma$  separation with N<sub>pe</sub>~10



## Beam test: Cherenkov angle resolution and number of photons



Beam test results with 2cm thick aerogel tiles: >4 $\sigma$  K/ $\pi$  separation





## How to increase the number of photons?



#### What is the optimal radiator thickness?



October 20, 2005

9th ICATPP Conference,





#### Beam tests



Photon detector: array of 16 H8500 PMTs

October 20, 2005

9<sup>th</sup> ICATPP Conf



#### Clear rings, little background







#### Focusing configuration – low momentum



Matching of indices: done for high momentum tracks (4GeV/c)
How is the overlapping of rings at lower momenta?



Good overlapping down to 0.6 GeV/c

9th ICATPP Conference, Como

Peter Križan, Ljubljana



#### Focusing configuration – momentum scan





Overlapp optimized at 4GeV/c-> OK at low momenta as well

October 20, 2005

9th ICATPP Conference, Como

Peter Križan, Ljubljana





# Optimisation of the dual radiator configuration









## Extension to many layers – expectations



1 -> 2 layers: a big jump in resolution per track

2 -> 3, 4: small steps

But: more radiator layers -> shallower optimal thickness minimum

-> of particular importance in the vicinity of threshold, where the number of photons is more important \_\_\_\_\_





October 20, 2005

9th ICATPP Conference, Como

Peter Križan, Ljubljana



## Aerogel production R&D





Going to thicker radiators: advantage only if aerogels have a high transmission length -> R+D Improved optical quality for n~1.05 hydrophobic aerogel with a new solvent and precursor



9<sup>th</sup> ICATPP Conference,



## Aerogel production R&D



Further optimization for n = 1.050 samples





## Aerogel production: multilayer samples





- 2 (or more) layers with different n
- layers attached directly at molecular level
- easy to handle
- Insensitive to possible surface effect







## Tiling of the radiator



Minimize photon yield losses at the aerogel tile boundary: hexagonal tiling scheme





- Cut into hexagonal shape from a square block
- Machining device: use "water-jet" thanks to hydrophobic nature



### Summary



- Proximity focusing RICH with aerogel as radiator is by now a well proven experimental technique.
- Increase the number of Cherenkov photons: employ radiators with multiple refractive indices. Idea successfully tested in beam tests.
- Aerogel production: transmission length improved, new cutting methods tested, multiple layer samples.
- R&D issues: development and testing of a multichannel photon detector for high mag. fields
- mass production of large aerogel tiles
- readout electronics







October 20, 2005



## Development and testing of photon detectors for 1.5 T



- Baseline: large area HPD of the proximity focusing type
- Backup: MCP-PMT





## HPD development



## 59mm x 59mm active area (65%), 12x12 channels





#### Ceramic HPD box

First tests carried out. Problems with sealing the tube at the windowceramic box interface.

Waiting for the next batch in September.



### Photon detector R&D – backup option: Burle MCP-PMT



#### BURLE 85011 MCP-PMT:

.multi-anode PMT with 2 MCPs
.25 μm pores
.bialkali photocathode
.gain ~ 0.6 x 10<sup>6</sup>
.collection efficiency ~ 60%
.box dimensions ~ 71mm square
.64(8x8) anode pads
.pitch ~ 6.45mm, gap ~ 0.5mm
.active area fraction ~ 52%





### Photon detector R&D – Burle MCP-PMT bench tests



Study uniformity of the sensitivity over the surface

count rates - all channels: charge sharing at pad boundaries

single channel response:uniform over pad areaextends beyond pad area (charge

sharing) October nference 







#### • BURLE MCP-PMT mounted together with an array of 12(6x2) Hamamatsu R5900-M16 PMTs at 30mm pitch (reference counter)



October 20, 2005

9<sup>th</sup> ICATPP Conference, Como





**Resolution and number of photons (clusters)** 

- $\sigma_9 \sim 13 \text{ mrad}$  (single cluster)
- number of clusters per track N ~ 4.5
- $\sigma_9 \sim 6 \text{ mrad (per track)}$
- -> ~ 4  $\sigma \pi/K$  separation at 4 GeV/c

#### **Open questions**

```
Operation in high magnetic field:
```

.the present tube with  $25\mu m$  pores only works up to 0.8T, for 1.5T need ~10 $\mu m$ Number of photons per ring: too small. Possible improvements:

- .increase active area fraction (bare tube 63%->85%)
- .increase the photo-electron collection efficiency

(from 60% at present up to 70%)

-> Extrapolation from the present data 4.5 ->8.5 hits per ring

 $\sigma_{g}$ : 6 mrad -> 4.5 mrad (per track)

-> >5  $\sigma \pi/K$  separation at 4 GeV/c

Aging of MCP-PMTs ?



### Read-out electronics: ASIC under development



Need high density front-end electronics. Need high gain with very low noise amplifiers. Deadtimeless readout scheme-> Pipeline.

#### Develop an ASIC for the front-end electronics

- Gain : 5 [V/pC]
- Shaping time : 0.15 [  $\mu$  s]
- S/N : 8 (@2000[e])
- Readout : pipeline with shift register
- Package : 18 channels/chip





Peter Križan, Ljubljana



### **Read-out electronics**



ASIC controled and read-out by a control board (for tests: can also be done with standard VME modules+level adapters)

- Detailed evaluation of the system is under way.
- Preparation of the read-out of an array of 3x3 HAPDs in a beam test.



Backup options use chips from Ideas: a VA/TA based system developed by the K2K group at KEK and VA64TAP+LS64







From typical values (single photon resolution 13mrad and 6 detected photons) we can estimate the Cherenkov resolution per track: 5.3mrad;

-> 4.3sigma p/K separation a 4GeV/c.

Illustration of PID performance: Cherenkov angle distribution for pions at 4GeV/c and 'kaons' (pions at 1.1GeV/c with the same Cherenkov angle as kaons at 4GeV/c). Details: NIM paper

Photon detector: array of 16 H8500 PMTs





#### How to design radiator tiles: check losses at the tile boundary.



- Scan with the beam across the tile boundary. As expected, the yield is affected over a few mm in the vicinity of the boundary.
- A simple model (all photons hitting the boundary get lost) accounts for most of the dependence

## Reduce the fraction of tracks close to tile boundaries and corners.